

2. Create Sync Site




This step is best carried out using the [Sync site wizard](#). These instructions are retained here for reference only.

The first step is to create the sync site that the store will be 'Active' on.


1. Got to Admin > Preferences



2. On the Synchronise tab on the left-hand side and click on the unlock button and enter the password



- 3. Add new Sync Site 
- 4. Enter **Sync ID**, **User name** and **Password**. Note that **URL** and **Hardware ID** are not required.



- 5. Click **OK**



- Keep track of the password! If it is lost you'll have to resynchronise from scratch!
- Sussol will typically have them securely stored if Sussol has done the sync setup.

Previous: **[1. Where to start](#)** | | Next: **[3. Create the Store](#)**

From:
<https://wiki.msupply.foundation/> - **mSupply Foundation Documentation**

Permanent link:
https://wiki.msupply.foundation/en:mobile:setup:server_side:create_sync_site

Last update: **2021/07/15 23:00**

