2024/05/21 18:22 1/3 2. Create Sync Site

2. Create Sync Site




This step is best carried out using the Sync site wizard. These instructions are retained here for reference only.

The first step is to create the sync site that the store will be 'Active' on.


1. Got to Admin > Preferences



2. On the Synchronise tab on the left-hand side and click on the unlock button and enter the password



- 3. Add new Sync Site
- 4. Enter Sync ID, User name and Password. Note that URL and Hardware ID are not required.



5. Click OK



- Keep track of the password! If it is lost you'll have to resynchronise from scratch!
- Sussol will typically have them securely stored if Sussol has done the sync setup.

2024/05/21 18:22 3/3 2. Create Sync Site

Previous: 1. Where to start | | Next: 3. Create the Store

From:

https://wiki.msupply.foundation/ - mSupply Foundation Documentation

Permanent link:

https://wiki.msupply.foundation/en:mobile:setup:server_side:create_sync_site

Last update: 2021/07/15 23:00

