

2. Create Sync Site


The first step is to create the sync site that the store will be 'Active' on.

 As of 2019-10-24, this process is manual. Each sync site needs to be created manually by an operator, one-by-one.

1. Got to Admin > Preferences




2. On the Synchronise tab on the left-hand side and click on the unlock button and enter the password



3. Add new Sync Site 

4. Enter **Sync ID**, **User name** and **Password**. Note that **URL** and **Hardware ID** are not required.



5. Click **OK**



- Keep track of the password! If it is lost you'll have to resynchronise from scratch!
- Sussol will typically have them securely stored if Sussol has done the sync setup.

Previous: [1. Where to start](#) | | Next: [3. Create the Store](#)

From:
<https://wiki.msupply.foundation/> - **mSupply Foundation Documentation**

Permanent link:
https://wiki.msupply.foundation/en:mobile:setup:server_side:create_sync_site?rev=1583105025

Last update: **2020/03/01 23:23**

