2025/12/15 04:45 1/3 2. Create Sync Site

## 2. Create Sync Site

The first step is to create the sync site that the store will be 'Active' on.




As of 2019-10-24, this process is manual. Each sync site needs to be created manually by an operator, one-by-one.


1. Got to Admin > Preferences



2. On the Synchronise tab on the left-hand side and click on the unlock button and enter the password



- 3. Add new Sync Site
- 4. Enter Sync ID, User name and Password. Note that URL and Hardware ID are not required.



## 5. Click OK



- Keep track of the password! If it is lost you'll have to resynchronise from scratch!
- Sussol will typically have them securely stored if Sussol has done the sync setup.

2025/12/15 04:45 3/3 2. Create Sync Site

Previous: 1. Where to start | | Next: 3. Create the Store

From:

https://wiki.msupply.foundation/ - mSupply Foundation Documentation

Permanent link:

https://wiki.msupply.foundation/en:mobile:setup:server\_side:create\_sync\_site?rev=1583105025

Last update: 2020/03/01 23:23

