User Guide - mSupply Mobile

11. Indicators

Indicators are used for reporting on the performance of drug regimens (programs). Each indicator value provides data on a metric for a particular store and regimen over a given period.

Indicators allow an organisation to analyse regimen performance and determine what further data should be collected, stored and processed.

11.1. Creating indicator data as a customer

If your store program has indicators associated with it you can add indicator data to the Supplier Requisition while creating it.

In the requisition window click the "Indicators" button in the top right corner. When data is added it will be sent to the supplier once the requisition is finalised.

Program: Order type:	HIV Nouveau type d 24/1/2020 mobile	Period: Supplier: Months Stock: Comment:	Avril - 2019 (04/09/19 - 05/16/19) Y				Items	Indicators
Entry Date: Entered By:			3			HIV		
Q Search	by item name or iter	n code Patients adultes	Patients enfants	Patient stable	Patient stable	Nouvelle	Nouvelle	Referes
e en en ip aren i		reçus	reçus	adulte	enfant	inclusion adulte	inclusion enfant	
	AZT/3TC/NVP	0	0	0	0	0	0	0
ière Ligne Régime ARV de	AZT/3TC/NVP AZT/3TC/EFV	0	0	0	0	0	0	0
Régime ARV de lère Ligne Régime ARV de lère Ligne Régime ARV de lère Ligne				-		-	-	

The indicators are configured on the server for each program, so what is seen will depend on what is has been centrally configured and what program you chose when creating the requisition.

11.2. Viewing indicator data as a supplier

When a corresponding Customer Requisition has been generated the indicator data can be viewed in the requisition by clicking the "Indicators" button in the top right corner.

Last update: 2020/01/30 en:mobile:user_guide:indicators https://wiki.msupply.foundation/en:mobile:user_guide:indicators?rev=1580428086

Note that in the Customer Requisition the indicator data can be viewed but not edited.

Previous: The Mobile Dashboard Next: Mobile Settings

From: https://wiki.msupply.foundation/ - mSupply Foundation Documentation

Permanent link: https://wiki.msupply.foundation/en:mobile:user_guide:indicators?rev=1580428086

Last update: 2020/01/30 23:48

